

# DICOM and NIfTI Data Standards in R

Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

Brandon Whitcher PhD CStat

Mango Solutions

London, United Kingdom

[www.mango-solutions.com](http://www.mango-solutions.com)

[bwhitcher@mango-solutions.com](mailto:bwhitcher@mango-solutions.com)

@Mangolmaging

12 June 2012 – useR!2012 Tutorial

- 1 Summary
- 2 Digital Imaging and Communications in Medicine
- 3 Neuroimaging Informatics Technology Initiative
- 4 Data Visualization
- 5 DICOM-to-NIfTI Conversion
- 6 Conclusions

- Content presented here provides users with the skills to manipulate DICOM / ANALYZE / NIfTI files in **R**.
  - Real-world data sets are used to illustrate the basic functionality of **oro.dicom** and **oro.nifti**.
  - S4 classes “nifti” and “anlz” enable further statistical analysis in **R** without losing contextual information from the original ANALYZE or NIfTI files.
  - Images in the metadata-rich DICOM format may be converted to NIfTI semi-automatically using as much information from the DICOM files as possible.
  - Basic visualization functions, are provided for “nifti” and “anlz” objects.
  - The **oro.nifti** package allows one to track every operation on a nifti object in an **XML**-based audit trail.

- There are several **R** packages that are able to access the DICOM / ANALYZE / NIfTI data formats:
  - **AnalyzeFMRI** (Bordier et al. 2009)
  - **fmri** (Polzehl and Tabelow 2007)
  - **Rniftilib** (Granert 2010)
  - **tracter.base** (Clayden 2010)

## Question #1

What are the (dis)advantages to having a single **R** package that performs input / output for medical imaging data?

## Question #2

Should **R** packages be discouraged from writing output in formats other than ANALYZE or NIfTI?

- The DICOM “standard” for data acquired using a clinical imaging device is very broad and complex.
  - Each DICOM-compliant file is a collection of fields organized into two two-byte sequences (group,element) that are represented as hexadecimal numbers and form a tag.
  - The (group,element) combination establishes what type of information is forthcoming in the file.
- There is no fixed number of bytes for a DICOM header.
- The final (group,element) tag should be the *PixelData* tag (7FE0,0010), such that all subsequent information is related to the image.

## DICOM

Digital Imaging and Communications in Medicine

<http://medical.nema.org>

Data element with explicit VR of OB, OF, OW, SQ, UT or UN:

```
+-----+
| 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 |
+-----+
|<Group-->|<Element>|<VR----->|<0x0000->|<Length----->|<Value->
```

Data element with explicit VR other than as shown above:

```
+-----+
| 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
+-----+
|<Group-->|<Element>|<VR----->|<Length->|<Value->
```

Data element with implicit VR:

```
+-----+
| 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
+-----+
|<Group-->|<Element>|<Length----->|<Value->
```

- Byte ordering for a single (group,element) tag in the DICOM standard.
- Explicit VRs store the VR as text characters in two bytes.

## readDICOMFile()

Accessing information stored in a single DICOM file

- The resulting object is a list with two elements: the DICOM header (hdr) and the DICOM image (img).
- Header information is organized in a data frame with six columns and an unknown number of rows.
- First five columns taken from DICOM header information (group, element, code, length and value) or inferred (name).
- (group,element) values are stored as character strings – not hexadecimal numbers.

## readDICOM()

Accessing multiple DICOM files in a single directory or directory tree

## create3D() and create4D()

Create arrays from DICOM headers / images

- Minimum input = “dcm” structure
- *PixelData* may be read on-the-fly
- Siemens MOSAIC format allowed
- 4D volumes may require additional information
  - nslices = ?
  - ntimes = ?



Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

- The ANALYZE format was originally developed in conjunction with an image processing system (of the same name) at the Mayo Foundation.
- A common version of the format, although not the most recent, is called ANALYZE 7.5.
- An ANALYZE 7.5 format image is comprised of two files, the “.hdr” and “.img” files, that contain information about the acquisition and the acquisition itself, respectively.

## The ANALYZE Format

.hdr = 348-byte binary file of header information

.img = binary flat file of images (multi-dimensional array)

- The NIfTI-1 data format is almost identical to the ANALYZE format, but offers a few improvements:
  - Merging of the header and image information into one file “.nii”
  - Re-organization of the 348-byte fixed header into more relevant categories
  - Possibility of extending the header information
- Discussions have begun on the NIfTI-2 data format.

## NIfTI

Neuroimaging Informatics Technology Initiative

<http://nifti.nimh.nih.gov>

Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

- “nifti”
  - Inherits from class “array”
  - Slots contain NIfTI header information
  - Basic methods: show(), nifti(), is.nifti(), as(<obj>, “nifti”)
  - Input / output: readNIfTI(), writeNIfTI()
  - Slot access: pixdim(), qform(), sform(), descrip(), aux.file(), audit.trail()
  - Additional classes: “niftiExtension” and “niftiExtensionSection”
- “anlz”
  - Inherits from class “array”
  - Slots contain ANALYZE header information

- **oro.nifti** offers three functions for visualization
  - `image()` = overloaded function for “anlz”, “array” and “nifti” classes
  - `overlay()` = extension of `image()` with `x` and `y` input parameters
  - `orthographic()` = mid-axial, mid-sagittal, mid-coronal views

## Interactive Visualization

### FSLView

<http://www.fmrib.ox.ac.uk/fsl/fslview/>

### MRICron

<http://www.cabiatl.com/mricro/mricron/>

### VolView

<http://www.kitware.com/products/volview.html>

Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

- **oro.dicom** and **oro.nifti** were designed to use as much information as possible from the metadata-rich DICOM format and apply that information in the construction of the NIfTI data volume.
  - Read in a single series using `dicomSeperate()`
  - `dicom2nifti()` converts the list of DICOM images into a multidimensional “nifti” object
  - `dicom2analyze()` converts the list of DICOM images into a multidimensional “anlz” object
  - Additional scripting in **R** is required to deal with multiple series

- Content presented here provides users with the skills to manipulate DICOM / ANALYZE / NIfTI files in **R**.
  - Real-world data sets are used to illustrate the basic functionality of **oro.dicom** and **oro.nifti**.
  - S4 classes “nifti” and “anlz” enable further statistical analysis in **R** without losing contextual information from the original ANALYZE or NIfTI files.
  - Images in the metadata-rich DICOM format may be converted to NIfTI semi-automatically using as much information from the DICOM files as possible.
  - Basic visualization functions, are provided for “nifti” and “anlz” objects.
  - The **oro.nifti** package allows one to track every operation on a nifti object in an **XML**-based audit trail.

Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

- <http://rigorousanalytics.blogspot.com>
- <https://r-forge.r-project.org/projects/rigorous>
- [rigorous-dicom@lists.r-forge.r-project.org](mailto:rigorous-dicom@lists.r-forge.r-project.org)
- [rigorous-nifti@lists.r-forge.r-project.org](mailto:rigorous-nifti@lists.r-forge.r-project.org)
- [rigorous-dcemris4@lists.r-forge.r-project.org](mailto:rigorous-dcemris4@lists.r-forge.r-project.org)

## Volume 44 in the *Journal of Statistical Software*

Special volume on “Magnetic Resonance Imaging in R”

- 13 articles on fMRI, DTI, DCE-MRI, etc.
- [www.jstatsoft.org/v44](http://www.jstatsoft.org/v44)

Summary

DICOM

Nifti

Visualization

Conversion

End

Bibliography

- C. Bordier, M. Dojat, and P. Lafaye de Micheaux.  
**AnalyzeFMRI**: an **R** package to perform statistical analysis on fmri datasets. 2009. URL <http://www.biostatisticien.eu/AnalyzeFMRI/>.  
Software: **R** Package, **AnalyzeFMRI**, version 1.1-12.
- J. Clayden. ***tractor.base**: A Package for Reading, Manipulating and Visualising Magnetic Resonance Images*, 2010. URL <http://CRAN.R-project.org/package=tractor.base>.  
**R** package version 1.5.0.
- O. Granert. ***Rniftilib**: R Interface to **NIFTICLIB** (V1.1.0)*, 2010. URL <http://CRAN.R-project.org/package=Rniftilib>. **R** package version 0.0-29.



Summary

DICOM

NIfTI

Visualization

Conversion

End

Bibliography

J. Polzehl and K. Tabelow. **fmri**: A package for analyzing fmri data. *RNews*, 7(2):13–17, 2007. URL [http://www.r-project.org/doc/Rnews/Rnews\\_2007-2.pdf](http://www.r-project.org/doc/Rnews/Rnews_2007-2.pdf).

- **Title:** Rigorous - DICOM Input / Output
- **Description:** Data input/output functions for data that conform to the Digital Imaging and Communications in Medicine (DICOM) standard, part of the Rigorous Analytics bundle.
- **Depends:** R ( $\geq 2.13.0$ ), utils
- **Suggests:** hwriter, oro.nifti ( $\geq 0.2.9$ )
- **License:** BSD
- **URL:** <http://rigorousanalytics.blogspot.com>

- **Title:** Rigorous - NIfTI+ANALYZE+AFNI Input / Output
- **Description:** Functions for the input/output and visualization of medical imaging data that follow either the ANALYZE, NIfTI or AFNI formats. This package is part of the Rigorous Analytics bundle.
- **Depends:** R ( $\geq 2.13.0$ ), bitops, graphics, grDevices, methods, utils
- **Suggests:** XML
- **Imports:** splines
- **Enhances:** fmri
- **License:** BSD
- **URL:** <http://rigorousanalytics.blogspot.com>